



Gas Chromatography- Mass Spectrometry

FM-GCMS-A100

Index

Sr. No	Title	Page no
1.	Introduction	2
2.	Features	2
3.	Specifications	3
4.	Applications	3
5.	Operations	4
6.	Accessories	17

1. Introduction

Gas Chromatography-Mass Spectrometry FM-GCMS-A100 supports an inlet temperature of up to 450°C, enabling a wide range of compound analysis. This is equipped with a turbo molecular pump, which maintains a stable vacuum environment essential for accurate analysis. Its column oven can achieve heating rates up to 120°C/min, facilitating faster analyses. Our GC-MS includes a dynode electron multiplier for superior sensitivity and detection of analytes.

2. Features

1. Versatile split injection modes
2. Integrated data processing system
3. Automated compound identification
4. High-speed scanning capability
5. Enhanced ion source efficiency

3. Specifications

Model	FM-GCMS-A100
Inlet Temperature	Max. 450°C
Pressure Range	0 to 100 psi (± 0.002 psi)
Pressure Control Mode	Electronic Pressure Control (EPC), supports Constant Volume (CV) and Constant Current (CC)
Split Mode	Split/Splitless, max. split ratio: 1000:1
Column Oven Temperature	Room temperature + 4°C to 450°C
Heating Rate	Up to 120°C/min
Temperature Programming	7 stages / 8 platforms
Auto Sampler	Optional
EI Source Ionization Energy	Adjustable from 5 eV to 250 eV
Mass Range	1.5 to 1000 amu
Resolution	Unit resolution (full width at half maximum)
Ion Source Temperature	100 to 350°C
Filament Emission Current	0 to 350 μ A
GC-MS Temperature Interface	Max. 450°C
Stability	± 0.10 amu over 48 hours
Sensitivity	Full scan, 1 pg OFN at m/z 272 with S/N $\geq 30:1$ (RMS)
Scan Rate	Up to 10,000 amu/s
Vacuum	Turbo molecular pump (67 L/s)
Detector	High-energy dynode electron multiplier
Power Consumption	3000W
Dimensions (L×H×W)	500 × 480 × 1200 mm
Weight	200kg

4. Applications

Gas Chromatography-Mass Spectrometry FM-GCMS-A100 is widely used in environmental monitoring, food safety, petrochemical analysis, and pharmaceutical research. It is essential for detecting pollutants, analyzing additives, and identifying trace compounds in various industries.

5. Operations

5.1 Analysis of phthalates

1) Analysis Condition

Column: DB-5MS, 30 m × 0.25 mm × 0.25 μm

Injection Mode: Splitless

Injection Time: 1.2 min

Carrier Gas Flow Rate: 1 mL/min (constant flow mode)

Inlet Temperature: 280°C

Ion Source Temperature: 250°C

Interface Temperature: 280°C

Solvent Rinse Time: 4 min

Oven Temperature: 60°C (hold for 4 min)

2) Phthalates standard spectra and parameters




Figure 1 Kinds of phthalates standard TIC

Phthalates component name, retention time, and characteristic ions

Ingredient Name	Retention time (min)	Target ion (m/z)	Reference ion (m/z)
Dimethyl phthalate (DMP)	7.83	163	77,133
Diethyl phthalate (DEP)	8.68	149	177,176
Phthalic acid diisobutyl (DIBP)	10.53	149	57,150
Dibutyl phthalate (DBP)	11.36	149	150,205

Phthalate, di(2 - methoxy) ester (DMEP)	11.68	59	149 and 104 haplotypes
Phthalate, di(4 - methoxy - 2 - pentyl) carbonate (BMPP)	12.40	149	85,167
Phthalate, bis (2-ethoxy)-ethyl ester (DEEP)	12.85	149	73,72
Dipentyl phthalate ester (DPP)	13.31	149	150,237
Phthalate, dihexyl ester (DHXP)	15.69	149	150,251
Phthalate, butyl benzyl phthalate (BBP)	15.83	149	91,206
Phthalic acid (2-butoxy)ethyl (DBEP)	17.40	149	57,193
Phthalate, dicyclohexyl phthalate (DCHP)	18.13	149	167,249
Phthalate, di (2 - ethyl)hexyl acrylate (DEHP)	18.35	149	167,279
Phthalate, di-n-octyl phthalate (DNOP)	21.15	149	150,279
O-dinonyl phthalate (DNP)	24.43	149	150,293

5.2 Analysis of polybrominated biphenyls, polybrominated diphenyl ethers

Analysis Condition

Column: DB-5MS, 15 m × 0.25 mm × 0.1 μm

Injection Mode: Splitless

Injection Time: 1.2 min

Carrier Gas Flow: 1.3 mL/min (constant flow mode)

Inlet Temperature: 280°C

Ion Source Temperature: 280°C

Interface Temperature: 310°C

Solvent Clearing Time: 3 min

Scan Range: 50–1000 m/z

Column Temperature Program:

- Initial Temperature: 100°C (hold for 2 min)
- Ramp Rate: 20°C/min to 240°C
- Ramp Rate: 30°C/min to 320°C (hold for 5 min)




Figure-2 Polybrominated diphenyl ethers standard solution TIC Figure (2mg / L)

Polybrominated diphenyl ethers component name, retention time and characteristic ions

Ingredient Name	Retention time/min	Reference ion
A polybrominated diphenyl ether monobromobiphenyl ether	0.373	250,248,141
Dibromobiphenyl ether dibromobiphenyl ether	2.392	328,326,168
The three polybrominated diphenyl ethers tribromobiphenyl ether	3.918	408,406,248
Tetrabromodiphenyl ether tetrabromobiphenyl ether	5.382	488,486,326
Pentabromodiphenylether pentabromobiphenyl ether	6.629	564,406,404
Six polybrominated diphenyl ethers hexabromobiphenyl ether	7.374	643,484,482
Seven polybrominated diphenyl ethers heptabromobiphenyl ether	8.251	722,562,456
Octa octabromobiphenyl ether	9.108	801,642,639
Nine polybrominated diphenyl ethers nonabromobiphenyl ether	9.947	881,721,719
Decabromodiphenyl ether decabromobiphenyl ether	10.752	959,799,797




Figure-3 part of the standard solution of polybrominated diphenyl TIC Figure (5mg / L)

Part of polybrominated biphenyls component name, retention time and characteristic ions

Ingredient Name	Retention time/min	Reference ion
Monobromobiphenyl monobromobiphenyl	0.472	234,232,152
Dibromobiphenyl dibromobiphenyl	1.356	312,310,152
Tribromobiphenyl tribromobiphenyl	1.844	392,390,230
Tetrabromodiphenyl tetrabromobiphenyl	2.503	470,310,308
Pentabromodiphenyl pentabromobiphenyl	2.927	550,390,388
Hexabromobiphenyl hexabromobiphenyl	3.456	628,468,466
Decabromodiphenyl decabromobiphenyl	10.303	944,783,781

5.3 Analysis of Pesticide Residues

1) Organophosphorus pesticide analysis

Column: DB-5MS, 30m x 0.25mm x 0.25µm

Injection port temperature: 230°C

Interface temperature: 260°C

Ion source temperature: 220°C

Oven temperature program:

- 50°C (hold for 2 min)
- Ramp at 25°C/min to 200°C (hold for 7 min)
- Ramp at 20°C/min to 280°C (hold for 5 min)

Carrier gas flow rate: 1 mL/min (splitless injection)

Injection time: 1 min

Scan range: 45–400 m/z

Solvent delay: 5 min




Figure-4 16 kinds of organic phosphorus standard solution TIC Figure

Dichlorvos:

- | | |
|------------------------|------------------------|
| 2. Diazinon | 11. A poisonous insect |
| 3. Fonofos | 12. Methidathion |
| 4. Chlorpyrifos-methyl | 13. B parathion |
| 5. Methyl parathion | 14. Carbophenothion |
| 6. Methyl Pirimiphos; | 15. Phosalone |
| 7. Fenitrothion | 16. Methyl Azinphos |
| 8. Malathion | |
| 9. Chlorpyrifos | |
| 10. Parathion | |

2) Organochlorine pesticide analysis

Of 13 organochlorine test conditions:

Column: DB-5MS, 30m x 0.25mm x 0.25µm

Injection port temperature: 260°C

Oven temperature program:

- Ramp from 90°C to 200°C at 20°C/min (hold for 10 min)
- Ramp from 200°C to 280°C at 35°C/min (hold for 3 min)

Carrier gas flow rate: 1 mL/min (splitless injection)

Injection time: 1 min

Interface temperature: 280°C

Ion source temperature: 240°C

Scan range: 45–400 m/z

Solvent delay: 4 min




Figure-5 13 Organochlorine standard solution TIC Figure

- | | |
|-----------------------------|---------------|
| 1. A-BHC | 9. 2,4'-DDD |
| 2. B-BHC | 10. Endrin |
| 3. Γ -BHC heptachlor | 11. 4'4'-DDD |
| 4,5. Aldrin. | 12. 2,4'-DDT |
| 6. Heptachlor EPOXIDE | 13. 4,4'- DDT |
| 7. 4,4'-DDD | |
| 8. Di reagents | |

5.4 Polycyclic aromatic hydrocarbons analysis

Analytical Conditions

Column: DB-5MS, 30m x 0.25mm x 0.25 μ m

Ion source temperature: 300°C

Inlet temperature: 280°C

Interface temperature: 280°C

Carrier gas flow rate: 1 mL/min (high purity He)

Injection method: Splitless

Injection time: 1 min

Oven temperature program:

- 50°C (hold for 1 min)
- Ramp at 25°C/min to 200°C
- Ramp at 8°C/min to 315°C (hold for 5 min)




Figure-6 16 PAHs the standard samples TIC Figure

5.5 The analysis of volatile organic compounds (VOC)

Headspace Conditions:

1) Pressure:

- Pressure: 0.12 MPa
- Inlet pressure: 0.07 MPa
- Pressure time: 0.3 minutes
- Injection time: 0.1 minutes
- Sample Chamber Temperature: 80 °C
- Valve Box Temperature: 80 °C
- Line Temperature: 90 °C

2) GC Conditions:

- **Column:** 60 m × 0.32 mm × 1.8 µm (HP-VOC)
- **Oven Temperature:**
40 °C (initial hold time: 2 minutes)
15 °C/min ramp to 100 °C (hold for 1 minute)
5 °C/min ramp to 180 °C
- **Column Flow:** 1.5 mL/min
- **Split Ratio:** 10:1
- **Inlet Temperature:** 180 °C
- **Carrier Gas:** High purity helium

3) MS Conditions:

- Interface Temperature: 180 °C
- Ion Source Temperature: 200 °C




Figure 7: The 14 VOC Standards TIC Figure

- | | |
|------------------------|-------------------------|
| 1) Ethanol | 8) n-Propyl acetate |
| 2) Isopropanol | 9) 4-Methyl-2-pentanone |
| 3) Butanone | 10) Toluene |
| 4) Ethyl acetate | 11) Acetic acid |
| 5) Isopropyl acetate | 12) n-Butyl acetate |
| 6) Benzene | 13) Ethylbenzene |
| 7) Glycol methyl ether | 14) Xylene |
| 15) Cyclohexyl ketone | |

5.6 Analysis of Formaldehydes

Analytical conditions:

1) Column:

- Column Type: DB-5MS
- Dimensions: 30 m × 0.25 mm × 0.25 μm

2) GC Inlet Conditions:

- Inlet Temperature: 25 °C
- Injection Volume: 1 μL
- Injection Mode: Splitless
- Splitless Injection Time: 1 minute
- Solvent Delay: 4 minutes
- Carrier Gas Flow Rate: 1 mL/min

3) MS Conditions:

- Interface Temperature: 280 °C
- Ion Source Temperature: 240 °C
- Scan Range: 50-300 amu

4) Oven Temperature Program:

- Initial Temperature: 80 °C (hold for 2 minutes)
- Ramp Rate: 30 °C/min to 250 °C (hold for 5 minutes)




Figure-8 formaldehyde derivatives TIC Figure

5.7 Azo dyes

Analysis Condition

1) Column:

- Column Type: DB-5MS
- Dimensions: 30 m × 0.25 mm × 0.25 μm

2) GC Inlet Conditions:

- Inlet Temperature: 260 °C
- Injection Mode: Splitless
- Injector Time: 1 minute
- Carrier Gas Flow Rate: 1.2 mL/min (High purity He)

3) MS Conditions:

- Ion Source Temperature: 200 °C
- Interface Temperature: 270 °C
- Scan Range: 35-350 amu

4) Oven Temperature Program:

- Initial Temperature: 60 °C (hold for 1 minute)
- Ramp Rate: 8 °C/min to 210 °C (hold for 10 minutes)
- Ramp Rate: 2 °C/min to 215 °C (hold for 1 minute)
- Ramp Rate: 20 °C/min to 260 °C (hold for 2 minutes)




Figure-9 23 kinds of azo dyes TIC figure

- | | |
|--|---|
| 1) 4-Amino-biphenyl | 14) 2-Methoxy-5-methylaniline |
| 2) Benzidine | 15) 4,4'-Alkylene-methyl-di-(2-chloroaniline) |
| 3) 4-Chloro-o-toluidine | 16) 4,4'-Diaminodiphenyl ether |
| 4) 2-Naphthylamine | 17) 4,4'-Diaminodiphenyl sulfide |
| 5) O-Aminoazobenzene | 18) o-Toluidine |
| 6) Toluene | 19) 2,4-Diaminotoluene |
| 7) Chloroaniline | 20) 2,4,5-Trimethylaniline |
| 8) 2,4-Amino-anisole | 21) Anisidine |
| 9) 4,4'-Diamino-diphenylmethane | 22) 2,4-Dimethylaniline |
| 10) 3,3'-Dichlorobenzidine | 23) 2,6-Xylene aniline |
| 11) 3,3'-Dimethoxybenzidine | 24) 4-Aminoazobenzene |
| 12) 3,3'-Dimethylbenzidine | |
| 13) 3,3'-Dimethyl-4,4'-diaminodiphenyl | |

5.8 The analysis of Tetrabromo bisphenol A, bisphenol A

Analysis conditions:

1) Column:

- Column Type: DB-5MS
- Dimensions: 30 m × 0.25 mm × 0.25 μm

2) GC Inlet Conditions:

- Inlet Temperature: 280 °C
- Injection Volume: 1 μL
- Injection Mode: Splitless
- Solvent Delay: 4 minutes
- Carrier Gas Flow Rate: 1.2 mL/min (High purity He)

3) MS Conditions:

- Interface Temperature: 300 °C
- Ion Source Temperature: 260 °C
- Scan Range: 60-650 amu

4) Oven Temperature Program:

- Initial Temperature: 100 °C (hold for 2 minutes)
- Ramp Rate: 30 °C/min to 310 °C (hold for 8 minutes)




Figure-10: Tetrabromobisphenol A, bisphenol A derivatives

5.9 Chlorophyll analysis

Analysis Condition

Column: DB-5MS, 30 m × 0.25 mm × 0.25 μm

Injection Type: Split injection, Split ratio: 10:1

Carrier Gas Flow Rate: 1 mL/min

Inlet Temperature: 240 °C

Interface Temperature: 270 °C

Ion Source Temperature: 240 °C

Scan Range: 50-320 amu

Oven Temperature Program: 100 °C (initial), ramp 10 °C/min to 240 °C (hold for 1 min).




Figure-11 Chlorophenol derivatives TIC

- 1) Pentachlorophenol
- 2) 2,3,4,5-Tetrachlorophenol
- 3) 2,3,5,6-Tetrachlorophenol
- 4) 2,4,6-Trichlorophenol

- 5) 2,3,5-Trichlorophenol
- 6) 2,4,5-Trichlorophenol

5.10 Analysis of dimethyl fumarate

Analysis Conditions

Column: DB-5MS, 30 m × 0.25 mm × 0.25 μm

Interface Temperature: 280 °C

Inlet Temperature: 280 °C

Ion Source Temperature: 230 °C

Carrier Gas Flow Rate: 1 mL/min

Injection Type: Splitless

Splitless Time: 1 min

Oven Temperature Program:

- 60 °C (initial)
- Ramp 5 °C/min to 100 °C
- Ramp 25 °C/min to 280 °C (hold for 10 min)




Figure-12: Dimethyl fumarate TIC Figure

5.11 Organotin analysis

Analysis Condition

Column: DB-5MS, 30 m × 0.25 mm × 0.25 μm

Ion Source Temperature: 250 °C

Inlet Temperature: 260 °C

Interface Temperature: 280 °C

Carrier Gas Flow Rate: 1.1 mL/min (high purity He)

Injection Mode: Splitless injection

Injection Time: 1 min

Solvent Delay: 2.5 min

Oven Temperature Program:

- 50 °C (hold for 1.5 min)
- Ramp 10 °C/min to 300 °C (hold for 4 min)




Figure-13: Kinds of organotin derivatives SIM scan spectra

- 1) Dimethyltin
- 2) Trimethyltin
- 3) Tributyltin
- 4) Dibutyltin
- 5) Phenyltin
- 6) Tributyltin
- 7) Diphenyltin

6. Accessories

Accessories	Quantity
Mechanical Pump	1 Pc
GC-MS Column	1 Pc
Network Cable	1 Pc
Series Port Cable	1 Pc
Power cable	5 Pcs
Corrugated pipe	1 Meter
Clamp hoop	2 Pcs
Seal ring	2 Pcs
Wrench - 5.5mm, 7mm	1 Pc
Wrench-10mm, 12mm	1 Pc
Wrench -8mm, 10mm	1 Pc
Tweezers	1 Pc
Cross Screwdriver - $\Phi 6 \times 120$, ph2	1 Pc
Cross Screwdriver- $\Phi 2 \times 50$, ph1	1 Pc
Straight Screwdriver	1 Pc
Monkey Wrench	1 Pc
Allen Wrench	1 Pc
Sleeve	1 Pc
15% Graphite /85% Vespel Ferrule	10 Pcs
Syringe	1 Pc
Column Cutter	1 Pc
Liner	1 Pc
Fuse	2 Pcs
Plastic Corrugated pipe	1.5 Meter
Sylphon Bellows Connector	1 Pc
Hose clamp	1 Pc
Ion Source Holder	1 Pc
Column measuring tool	1 Pc
Stainless Steel Tube	7 Meter
1/8" Ferrule	1 Pc
Valve Connector	2 Pcs
O-ring - 2.5 $\times$ 1.8	5 Pcs
O-ring - 1.8 $\times$ 1.8	1 Pc
O-ring	5 Pcs
Glass Wool	1 Pc
Front Ferrule	3 Pcs
Back Ferrule	3 Pcs
Injection port Wrench	1 Pc
Glass Liner	1 Pc
Moisture Trap	1 Pc
Oxygen Trap	1 Pc
Trap Bracket	1 Pc

Ribbon	8 Pcs
Injection Septum	1 Pc
Ethernet switch	1 Pc



Fison Instruments Ltd 272 Bath Street Glasgow G2 4JR UK
Email: info@fison.com | Website: www.fison.com